




WILD TYPE SARS-COV-2 SURROGATE VIRUS CE NEUTRALISATION ASSAY

Benefits

- » Qualitative detection of neutralising antibodies against SARS-CoV-2 Wild Type variant
- » Enables identification of a protective immune response to SARS-CoV-2 Wild Type
- » Aids in the stratification of at risk groups where additional booster vaccinations may be required
- » Accurate assessment of herd immunity and measurement of vaccine efficacy
- » Capable of detecting neutralising antibodies in an isotype and species independent manner
- » Rapid results (1.5 hours) without the need for Biosafety level 3 requirements
- » Specificity of 100% and sensitivity of 98% compared to conventional virus neutralisation test
- » Direct sample addition to Biochip well. No off-board sample incubation required

Available on the Evidence Investigator, the Randox CE marked Wild Type SARS-CoV-2 Surrogate Virus Neutralisation assay utilises patented Biochip Array Technology (BAT) to detect neutralising antibodies which inhibit the interaction between SARS-CoV-2 Spike Receptor Binding Domain (RBD) and the cell surface Angiotensin-converting enzyme 2 (ACE2).



www.randox.com

marketing@randox.com

y D in CE


Copyright © 2021 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151682708. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representativew





EVIDENCE INVESTIGATOR Semi-Automated Analyser

The Evidence Investigator is a compact semi-automated benchtop analyser ideal for medium throughput laboratories.



ORDERING INFORMATION	
Wild Type SARS-CoV-2 Surrogate Virus Neutralisation assay	EV4454
Evidence Investigator Analyser	EV3602

www.randox.com

marketing@randox.com



Copyright © 2021 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151682708. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representativew