

HbA1c

The RX series world leading test menu welcomes the addition of direct HbA1c testing on the RX modena, RX imola and RX daytona+ chemistry analysers. This new immunoturbidimetric HbA1c test enables laboratories to expand their testing capabilities onto one single platform; offering cost savings through consolidation.

KEY BENEFITS OF DIRECT HbA1c TESTING ON THE RX SERIES



Single assay testing allows faster recovery times



No offline preparation stage ensures accurate sample & QC recovery



No incubation step allows samples to be run immediately



No offline calculations ensures faster recovery times

RX SERIES CLINICAL CHEMISTRY ANALYSERS



RX daytona+

- 450 tests per hour including ISE
- Low water consumption of 5 litres per hour
- STAT functionality
- Dual 5 speed stirrers
- Automatic wake and sleep mode functionality



RX imola

- 560 tests per hour including ISE
- Low water consumption of 18 litres per hour
- STAT functionality
- Liquid level sensor, crash bubble and crash detection
- Automatic wake and sleep mode functionality



RX modena

- 1200 tests per hour including ISE
- Interactive touchscreen technology
- Automatic standby and autostart functionality
- Economic water consumption of less than 20 litres per hour
- Continuous loading hatch for STAT samples

RANDOX

+44 (0) 28 9442 2413 • therxseries@randox.com • randox.com/hbalec-rxseries
Randox Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, United Kingdom, BT29 4QY

Product availability may vary from country to country. Please contact your local Randox representative for information.