



Anti- L-Selectin Recombinant Antibody Fragments

Target:	Anti- L-Selectin
Clone name:	MC39A10
Catalogue Number:	RAF9258
Description:	Recombinant single-chain variable fragment (scFv) ¹ obtained from Sheep and expressed in an E. Coli to bind against L-Selectin antigen.
Activity/ Specificity:	Tested positive against L-Selectin antigen. Cross-reactivity checked against a panel of known cross-reactants and non-specific antigens.
Applications:	These fragments contain His and c-Myc fusion tags which may be used for detection or immobilisation. ² Recombinant antibody fragments are suitable for use in ELISA immunoassays, biosensor applications, western blots, immunohistochemistry, flow cytometry, immunoaffinity purification and most other immunological methods*.
Size:	Approximately 35 kDa
Quantity:	1 mg
Concentration:	Typically >1mg/ml
Purity:	>90% assessed by SDS-PAGE.
Storage:	These fragments are stable at 4°C. It is recommended that for storage over extended periods they are kept at –20°C and should not be subject to repeated freeze-thaw cycles.
Buffer:	1x PBS containing 0.09% sodium azide preservative.
Dilution Factor:	To be determined by end-user.

* This product as supplied is intended for research applications only. It is not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox Life Sciences.

References:

- Weisser NE, Hall JC. (2009) Applications of single-chain variable fragment antibodies in therapeutics and diagnostics. Biotechnol Adv. 27(4): 502-520.
- 2. Terpe, K. (2003) Overview of tag protein fusions: from molecular and biochemical fundamentals to commercial systems. *Applied Microbiol Biotechnol.* 60(5):523-33.

Randox Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom T +44 (0) 28 9442 2413 F +44 (0) 28 9445 2912 E marketing@randox.com ₩ randox.com



Randox Biosciences is part of Randox laboratories W randoxbiosciences.com

Directors: Dr S P FitzGerald CBE FREng BSc PhD DSc (Hons) FIAE Man. R P Kelly BSc FCA.